



# TRANSIL Brain Absorption:

A novel in vitro model for predicting CNS penetration

Dr. Hinnerk Boriss

High Throughput Screening for Biological Systems  
March 16<sup>th</sup>, 2009

# Outline

- Significance of the Blood-Brain Barrier
- Key criteria for compound optimization
- What can we measure with existing models?
- Advantage of estimating both  $f_u$  brain and logBB

# TRANSIL Brain Absorption


## Key Features

- Prediction of brain availability ( $\log BB$  and  $f_u$  brain)
- Prediction of concentration at receptor site
- Classification of CNS+ /CNS- compounds
- Fast: 2 min incubation time
- Easy: First ready-to-use assay kit


# Significance of Blood-To-Brain Absorption

- Blood Brain Barrier protects the brain from xenobiotics while maintaining metabolic functions
- Significant hurdle for CNS targeted drugs
- Important protection from drug side-effects
- CNS-accessible chemical space much smaller in comparison to other organs

# Brain Barriers




# Neurovascular Unit



- Endothelial Cells
  - Interface btw blood and brain
  - Single cell encloses lumen
  - Tight junction holds cell together
- Tight junctions
  - Maintain structure
  - Signal transduction
  - Regulate paracellular diffusion
- Pericyte Cells
  - Flat contractile cells
  - Regulate endothelial cell proliferation
  - Induces BBB tightness
- Extracellular Matrix
  - 20% of brain
  - Organized mesh of secreted proteins
- Astrocytes
  - Surround 99% of BBB
  - Influence endothelial structure
  - Induce BBB tightness
- Interstitial Fluid
  - Receptor biophase

# Brain Compartments & Methods



# Key Issues

- Key question: will sufficient drug reach the site of action?
- Key issue: ISF if the receptor biophase

$$\text{ISF} \approx C_u\text{brain} = C_t\text{brain} \times f_u\text{brain}$$

## Key Criteria for Brain Absorption

- CNS targets
  - What is the concentration at the receptor site?
- CNS side effects
  - How much compound gets into brain?
  - How much compound is freely available to cause havoc in brain?

# What Can We Measure?

	Conc.	$f_u$
• Brain/Plasma distribution (logBB)		
– <i>in vivo</i>	+	-
– IAM	+	-
• Uptake rates		
– PAMPA	0	-
– MDCK, caco-2, BCEC, etc	0	-
• PS-product		
– <i>in situ</i> perfusion	0	-
• Brain free fraction		
– <i>in situ</i> dialysis	-	+
– Microdialysis against brain tissue	-	+

# Model Limitations

- logBB confounded by
  - Non-specific tissue binding
  - Efflux and influx transport
  - Brain metabolism
- Rate measurements (logPS, PAMPA, MDCK, etc) yield little information on brain conc.
- Plasma/CSF and Brain/SCF barrier different from Plasma/Brain barrier
- *in vivo*  $f_u$  determination expensive and laborious

# Classification Schemes


<b>Assay</b>	<b>CNS+</b>
MDR1-MDCKII	$P_{app} > 14 \cdot 10^{-6} \text{ cm/s}$
PAMPA-BBB	$P_e > 4 \cdot 10^{-6} \text{ cm/s}$
<i>in vivo</i> B/P ratio	$B/P > 0.3$
<i>in vivo</i> perfusion	$\text{Log PS} > -2.5$
CNS rules	$\text{PSA} < 60 - 70 \text{ \AA}^2$ $1 < \log D_{7.4} < 3$ $\text{MW} < 450$

## Free Drug Hypothesis and logBB

	A	B	C	D
$C_t$ blood	100	100	100	100
Blood $f_u$ [%]	1	10	1	10
Brain $f_u$ [%]	1	1	10	10

$C_u$ brain [ng/ml]	1	10	1	10
$C_t$ brain [ng/g]	100	1000	10	100
B/P ratio	1	10	0.1	1


# Brain Free Fraction and logBB



# BBB Models




Assay	Permeability		Equil. Dist.	$f_u$
	passive	active		
• In silico	☹️	-	☺️	-
• MDCK cells	☹️	☺️	-	-
• BCEC/astrocyte co-culture	☺️	☺️	-	-
• PAMPA	☺️	-	-	-
• TRANSIL Brain Absorption	☺️	-	☺️	☺️

# TRANSIL Brain Absorption I





# TRANSIL Brain Absorption II

- Porcine brain lipids
- Reconstituted membrane vesicles on 10  $\mu\text{m}$  beads
- Titration of membrane affinity




# Membrane – Small Molecule Interaction




# TRANSIL Brain Absorption A Hybrid Model

## Prediction based on:

- Brain membrane affinity }  $f_u$
  - Polar surface area
  - .....
  - Protein binding
- permeability
- equilibrium distribution
- 
- A diagram illustrating the relationship between input factors and equilibrium distribution. It features a list of three factors: "Brain membrane affinity", "Polar surface area", and "Protein binding". The first two factors are grouped by a right-facing curly bracket labeled
- $f_u$
- . A dotted line is placed between "Polar surface area" and "Protein binding". A larger right-facing curly bracket groups the
- $f_u$
- group and the "Polar surface area" factor, with the label "permeability" positioned to its right. A final, even larger right-facing curly bracket groups the
- $f_u$
- group, "permeability", and "Protein binding", with the label "equilibrium distribution" positioned to its right.

# Validation Data


## Compound Composition



Atenolol  
 Chlorambucil  
 Cimetidine  
 Clonidine  
 Ibuprofen  
 Indomethacin  
 Lidocain  
 Nalidixic acid  
 Phenylbutazone  
 Pindolol  
 Procaine  
 Propranolol  
 Pyrilamine  
 Quinidine  
 Ranitidine  
 Salicylic acid  
 Theophylline  
 Tolbutamide


Acetaminophen  
 Acetylsalic acid  
 Amitryptilline  
 Caffein  
 Carbamazepine  
 Chlorpromazine  
 Desipramine  
 Diphenhydramie  
 Eserine  
 Fluoxetine  
 Imipramine  
 Mianserin  
 Phenazone  
 Phenytoin  
 Promazine  
 Propofol  
 Thioperamide  
 Thioridazine  
 Trifluoperazine  
 Verapamil

# Validation Data: in vivo logBB



# Linear Model of Validation Data



Validation Data



Variable	<i>p</i> -value
$\log MA_{\text{brain}}$	$6.9 \times 10^{-7}$
$\log K_{B/F}$	$2.1 \times 10^{-4}$
PSA	$1.2 \times 10^{-8}$


$$\log BB = a \cdot \log MA_{\text{brain}} - b \cdot \text{PSA} - d \cdot \log K^{B/F} + \text{const}$$

# Key logBB Predictors: $\log\text{MA}_{\text{brain}}$ and PSA




# Leave One Out Cross Validation


## Bland-Altman Plot of Prediction Error



## Prediction of Brain Free Fraction



# CNS Classification



# Applications

- Primary screening in CNS projects
  - Filter: reduce library size by 30 to 60%
  - Build: support CNS specific library design
  - Decrease false negative rate
- Lead optimization in CNS projects
  - Optimize brain penetration
  - Cut cost by more than 50%
- Lead optimization in non-CNS projects
  - Optimize toxicology profile



# TRANSIL Brain Absorption Key Benefits

- Prediction of brain availability ( $\log BB$  and  $f_u$  brain)
- Classification of CNS+ /CNS- compounds
- Screening cost reduction by at least 50%
- Fast: 2 min incubation time
- Easy: First ready-to-use assay kit

# Contact Details

## **Sovicell GmbH**

Deutscher Platz 5b

04103 Leipzig

Germany

Tel: +49-341-52044-0

Fax: +49-341-52044-12

[info@sovicell.com](mailto:info@sovicell.com)

[www.sovicell.com](http://www.sovicell.com)

